

CHLORENOL LEAFLET FOR THE USER

Appearance of the preparation : liquid.

Instructions for use :

Working solutions and exposure times:

This preparation is a concentrate. It is used to prepare working aqueous solutions. It is prepared by soaking prior to disinfection of surgical instruments, laboratory glassware and patient-care equipment.

The disinfection of instruments is performed in a lidded tub. The solution should cover the instruments placed in the tub, and fill their inside openings.

The medical instruments should be thoroughly cleaned (from blood, tissues, proteins) prior to disinfection.

Following disinfection, the instruments should be rinsed thoroughly with clean or sterile (if required) water.

The disinfection solution may be used once within the 8-hour work day. The container of solution should always be covered with a lid to prevent evaporation. The instruments, prior to being placed in the solution, should be not only cleaned but also dry so as to prevent dilution.

The disinfection of laboratory equipment should also be performed in lidded tubs. For the patient-care equipment buckets and other containers may be used, also lidded.

Upon expiration of the contact time, both laboratory and patient-care equipment should be rinsed with potable water.

Disinfection should be performed in ventilated/well ventilated rooms.

Working solutions and exposure times:

- 3% working solution for 60 minutes (bactericidal and fungicidal effects),
- 5% working solutions for 30 minutes (bactericidal and fungicidal effects),
- 10% working solution for 15 minutes (bactericidal and fungicidal effects).

Number and duration of treatments : depend on the frequency of manipulations and the adopted disinfection practice.

First aid measures:

Treatment is symptomatic.

In case of inspiration: Ensure access to clean air. In case of discomfort, seek medical aid.

In case of contact with skin: Remove contaminated clothing immediately and wash it prior to the next use. Rinse any contaminated skin areas with soap and water. In case of irritation, seek medical assistance.

In case of contact with eyes: Rinse immediately under running water for 10-20 minutes with wide open eyelids so as to ensure the water getting under the eyelids. Medical advice should be sought.

In case of ingestion: Immediately rinse the mouth with water without ingesting it. Do not induce vomiting. Immediately seek medical aid and show the physician the packaging or the label of the packaging.

During use :

Personal protective equipment : safety glasses, gloves and protective clothing.

Gloves : latex or neoprene.

Avoid any contact with the eyes and skin. Do not breathe vapour.

Transportation :

This preparation can be transported by all types of transport in compliance with the rules for transporting goods appropriate for these types of transport. Do not transport along with food and drinks.

Storage :

Store in places protected from direct sunlight, away from heaters and open fire, at 5 to 30 .
Store in the original packaging, tightly sealed, away from the reach and sight of children.

Decontamination : In case of an incident and/or accidental spillage, take appropriate measures to localize and limit it, storing temporarily the accumulated amount of the preparation in special, tightly sealed and marked containers; afterwards, deliver it to the holders of license under art.37 of Waste Management Act (SG No 86/2003). Undertake measures for not allowing surface and ground water contamination, soil, or discharge into the sewage.

Appropriate firefighting agents: dry powders, dispersed stream of water, 2, alcoholic resistant foam.

Volumes: 200 ml

Batch No: see the pack

Date of manufacture: see the pack

Date of expiry: see the pack